

## Insert pro Pandemic

Původní návod k obsluze

### Technický popis

Výrobek je určen k třídění a skladování komponent deskové hry **Pandemic** v původní krabici základní hry včetně rozšíření Nové hrozby, Laboratoř a Stav ohrožení.

### Popis jednotlivých částí výrobku

Celá sestava se skládá z několika menších krabiček a pořadačů, které slouží k uskladnění komponent deskové hry. Popis uložení komponent je v **Appendixu 1**.

- 1 x Hlavní uložení pro přípravné karty
- 1 x Hlavní uložení pro hrací karty a žeton úspěchů
- 1 x Hlavní uložení Petriho misek, lahviček, kostky, figurek, atd.
- 1 x Krabička pro centra WHO a jiné budovy
- 1 x Krabička na žetony ze Stavů ohrožení

### Pokyny pro využívání výrobku

Výrobek je určen výhradně ke skladování a využívání komponent určené deskové hry. Výrobek není určen k dlouhodobému vystavování UV záření, teplotám nad 60 °C, mechanickému namáhání a rázovým silám.

### Pokyny pro údržbu výrobku

Z výrobku vyjměte všechny komponenty deskové hry a otřete jej vodou navlhčeným hadříkem.

### Pokyny k recyklaci

Veškerý materiál lze i po ukončení smyslu svého účelu recyklovat, a to dle následujících instrukcí.



Krabice: PĀP



Obalové fólie: PE-LD



Komponenty insertu: PET



Dokumentace: PĀP

**Prohlášení o shodě** je ke stažení včetně instrukcí na: <https://www.kalkared.eu/navody-a-dokumentace/>

## Insert for Pandemic

Original operational manual

### **Technical description**

Product is made to separate and store parts of a board game **Pandemic** in an original box of a core game with expansions On the Brink, Laboratory and State of the Emergency.

### **Product components description**

Assembly is made of several smaller boxes and trays which serves to store parts of a board game. How to store parts is described in **Appendix 1**.

- 1 x Main tray for setup cards
- 1 x Main tray for game cards and achievement tokens
- 1 x Main box for petri dishes, bottles, dice, pawns, etc.
- 1 x Box for WHO center and other buildings
- 1 x Box for State of Emergency tokens

### **How to properly use this product**

Product is exclusively made to store and use parts of designated board game.  
Product is not suitable to withstand long term UV radiation, temperatures above 60 °C, mechanical stress and shock loads.

### **Maintenance instructions**

Remove all parts of the board game from product and wipe it with moistened rag.

### **Recycle instructions**

All material can be recycled after its purpose comes to an end by following instructions.



Box: PAP



Packing foil: PE-LD



Insert components: PET



Documentation: PAP

**Declaration of conformity** can be downloaded with instructions at: <https://www.kalkared.eu/navody-a-dokumentace/>